



Optical detection of KCl vapor and atomic K releases from biomass fuels combusted in single particle reactor

Tapio Sorvajärvi and Juha Toivonen

Tampere University of Technology
Department of Physics, Optics Laboratory


Background

- KCl accelerates high temperature corrosion and reduces the heat transfer between flue gas and superheated water vapor
- KCl reduction techniques require real-time monitoring
- The study of KCl formation kinetics needs sensitive diagnostics




Optical Spectroscopy

- The study of the interaction between radiation and matter as a function of wavelength.





- Benefits in metrology:
 - Sensitivity
 - Wide dynamic range
 - Molecule specific
 - Possibility of remote and *in-situ* sensing


Absorption of KCl and K



Potassium chloride (KCl)



Potassium (K)




CPFAAS technique

(Collinear Photofragmentation and Atomic Absorption Spectroscopy)


Two lasers

UV laser emitting 1 ns pulses



Potassium chloride (KCl)

IR laser emitting cw-light




Potassium (K)




CPFAAS technique


1. IR laser transmission




2. KCl fragmentation



3. Vapor relaxation



CPFAAS technique



$$\frac{I_k}{I_{k0}} = \exp \left[- \left(1 - \exp \left(-x_{KCl} \frac{p}{k_b T} \sigma_{KCl} L \right) \right) \frac{E_p \sigma_K}{hc / \lambda_p A_p} \right]$$

$$x'_{KCl} = -\ln \left[1 + \ln \left(\frac{I_k}{I_{k0}} \right) \frac{A_p}{E_0} \frac{hc / \lambda_p}{\sigma_K} \frac{M_{KCl}}{N_a \sigma_{KCl}} \frac{1}{L} \right] \quad [x'_{KCl}] = \frac{g}{m^3}$$

Tunable parameters: —————

Physical constants: —————

Sample specific quantity: —————





Applications

- Determination of vapor pressures in thermodynamic equilibrium
- Real-time analysis of KCl release from combustion gas and flame
- Monitoring of KCl concentration in flue gases





Experiments - Single particle reactor

[4]





Results

20 mg of bagasse
950 C 3 % O₂





Results - Bagasse

950 C 20 mg





950 C 45 mg





Results - Bagasse



800 C 20 mg



950 C 20 mg




1100 C 10 mg




Results – Different Fuels


Bagasse 950 C 20 mg



Reed 950 C 22 mg




Straw 950 C 22 mg



Results – K release

Reed 950 C & 3 % O₂



Conclusions

- CPFAAS technique was applied to measure KCl and K releases from combusted biomass fuels
- KCl release was different from different fuels
- KCl release respect to CO₂ was found to depend on the combustion temperature
- K release profile changed as the mass of the sample increased.



Acknowledgements

The authors are grateful for the support by Doctoral Programme of TUT's President. The authors acknowledges the Process Chemistry Centre at Åbo Akademi University for providing facilities to the single particle reactor measurements and Aerosol Physics Laboratory at Tampere University of Technology for providing CO₂ data from single particle reactor experiments.

This work has been partly carried out within the consortium FUSEC (2011 – 2014) with support from National Technology Agency of Finland (Tekes), Andritz Oy, Metso Power Oy, Foster Wheeler Energia Oy, UPM-Kymmene Oyj, Clyde Bergemann GmbH, International Paper Inc. and Top Analytica Oy Ab.

